

High-end fluorescent microscope with ultrafast
precise liquid autofocus and cell counting software

MUCH MORE THAN A SIMPLE SOMATIC CELL COUNTER

Lactoscan Somatic Cell Counter

Embedded printer

10" fully functional
touch screen tablet

Test results within 20 to 60 seconds



SPECIFICATIONS

Analysis
time: < 1 minute

Capacity: > 60 tests per hour


Measuring
range: 10 000 to 10 000 000 cell/mL

Deviation: SCC (cell/mL)	CV%
100 000	5%
400 000	3%
600 000	2%

New

**DETECTION OF
YEAST MASTITIS!**

Accurate Somatic Cell Counting in three simple steps:



KEY FEATURES:

- Simple in operation
- Fast somatic cell analysis
- Accurate and reliable results
- Portable and compact design
- Low maintenance cost
- One year full warranty

Results are color coded



- ! Clinical Mastitis > 400 000 cell/mL
- ! Subclinical Mastitis 200 000 ÷ 400 000 cell/mL
- Healthy < 200 000 cell/mL
- ! Yeast mastitis

LACTOSCAN SCC STARTER KIT



**THE WORLD'S
LOWEST COST
PER TEST OF IMAGE
FLUORESCENCE
SOMATIC
CELLS COUNT**

This project is
co-funded by the
European Union



Improve your herd's
udder health

See with your eyes
what the others don't see

ENVIRONMENTAL CONDITIONS:

Ambient temperature: 10 - 40 °C
Milk temperature: 1- 40 °C
Relative humidity: 30 - 80 %

ELECTRICAL PARAMETERS:

AC voltage: 220V / 110V
DC voltage: 12V
Power Consumption:
30W max

MECHANICAL PARAMETERS:

Dimensions
(W x D x H):
390 x 300 x 260 mm
Weight < 6 kg

Plastic
cover box

buy online on
www.lactoscan.com




Detects subclinical
and clinical mastitis in
single teat of the udder



The iceberg of mastitis


**CLINICAL
MASTITIS**

**SUBCLINICAL
MASTITIS**




CCD camera Ultrafast precise liquid lens autofocus Fluorescent microscope


Somatic
cell images



3 working modes




Indication of different mastitis and pathogens




Normal milk 100x10³ cell/mL Staphylococcus mastitis Streptococcus mastitis Yeast mastitis Mycoplasma mastitis

The Principle of LACTOSCAN SCC



Unique, 3D, multi-image,
patent application protected,
sequential scanning process




LACTOSCAN SCC

compare to




DMSCC

Direct fluorescent,
low magnification
microscopic
somatic cell counter



Image

Direct
Microscopic
Somatic Cell
Counting



	LACTOSCAN SCC	DMSCC
SCC (cell/mL)	*CV %	*CV %
100 000	5 %	7 %
400 000	3 %	5 %
600 000	2 %	4 %

* Coefficient of Variation

Comparison of deviation

MILKOTRONIC LTD

Bulgaria, 6000 Stara Zagora
167, Tsar Simeon Veliki, A, 25
Tel: + 359 457 670 82,
Fax: + 359 457 670 83

e-mail: salesdpt@lactoscan.com
Eng, Spa, Port: +359 886 841 389
e-mail: sales2@lactoscan.com
Eng, Rus: +359 886 841 395
e-mail: mainservice@lactoscan.com
Eng, Rus: +359 888 811 017

e-mail: m.mihneva@lactoscan.com
Eng: + 359 885 597 232
e-mail: milkanalyzers@lactoscan.com
Eng, Rus: +359 896 753 803

distributors: